# EchelonVT

ECVT22L11S to ECVT30R84S  
ECVT22L11SD to ECVT30R84SD

EN	Instructions for Use	1
FR	Instructions d'utilisation	17
DE	Gebrauchsanweisung	33
IT	Istruzioni per l'uso	49
ES	Instrucciones de uso	65

endolite

# 1 Description and Purpose

## Application:

These instructions are for use by the practitioner.

EchelonVT is to be used exclusively as part of a lower limb prosthesis.

This device provides limited self-alignment of the prosthesis on varied terrain and following footwear changes. It is intended to improve on postural sway and symmetry while easing abnormal pressures at the socket interface. The device also includes a spring element capable of absorbing axial and rotational shocks intended to reduce shear stresses at the residuum/socket interface. This device is recommended for users that have the potential to achieve Activity Level 3 that may benefit from enhanced stability and an increase in confidence on uneven surfaces. The additional shock absorption and energy return supplied by the axial spring aims to improve comfort and performance.

There are exceptions and in our recommendation we want to allow for unique, individual circumstances. There may also be a number of users of Activity Levels 2 and 4\* who would benefit from the enhanced stability offered by EchelonVT, but this decision should be made with sound and thorough justification.

(\*maximum user weight 100kg and always use one higher spring rate category than shown in the Spring Set Selection table.)

- ⚠ **To minimise the risk of slipping and tripping, appropriate footwear that fits securely onto the footshell must be used at all times.**
- ⚠ **Be aware of finger trap hazard at all times.**

## Contra-indications:

This device may not be suitable for level 1 individuals or for competitive sports events, as these types of users will be better served by a specially designed prosthesis optimized for their needs.

Intended for single user.

Ensure that the user has understood all instructions for use, drawing particular attention to the section regarding Maintenance.

## Spring Set Selection

		User Weight								kg (lbs)
		44-52 (100-115)	53-59 (116-130)	60-68 (131-150)	69-77 (151-170)	78-88 (171-195)	89-100 (196-220)	101-116 (221-255)	117-125 (256-275)	
Activity		1	2	3	4	5	6	7	8	Foot Spring set
	3		•	••	•••	••••				
		1	2	3	4					

## Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

### Note:

If in doubt choosing between two categories, choose the higher rate spring set.

Foot Spring set recommendations shown are for transtibial users.

For transfemoral users we suggest selecting a spring set one category lower, refer to fitting advice Section 7 to ensure satisfactory function and range of movement

**Order example:**

<b>ECVT</b>	<b>25</b>	<b>L</b>	<b>N</b>	<b>3</b>	<b>2</b>	<b>S</b>
	Size	Side (L/R)	Width* (N/W)	Spring Set Category	Axial Spring	Sandal Toe

Available from size 22 to size 30:  
 ECVT22L11S to ECVT30R84S  
 ECVT22L11SD to ECVT30R84SD  
 (add 'D' for a dark tone foot shell)



\*Sizes 25-27 only. For all other sizes, omit the Width field.

e.g. ECVT25LN12S, ECVT22R34S, ECVT27RW33SD

## 2 Construction

**Principal parts:**

- Hydraulic Body and Pylon Assembly (Aluminium/St. Stl./Titanium)
- Carrier Assembly (Aluminium/St. Stl.)
- Heel & Toe Springs (e-Carbon)
- Spring attachment screws (St. Stl./Titanium)
- Glide Sock (UHM PE)
- Foot shell (PU)

---

## 3 Function

EchelonVT comprises a hydraulic body assembly containing adjustable hydraulic valves. The valves can be independently adjusted to increase and reduce hydraulic resistance to plantarflexion and dorsiflexion.

The proximal part of the hydraulic body assembly contains a spring component that allows the pyramid tube to deflect and/or rotate through a limited range under an axial and/or torsional load. When the load is removed the spring returns to its original position. The effect of this movement is to attenuate any shock forces that would otherwise be transmitted to the residuum/socket interface both in axial and rotational directions.

The hydraulic body assembly is connected to a carrier assembly via two pivot pins. Heel and toe springs are attached to the carrier using titanium and stainless steel screws. The foot is wrapped in a UHM PE sock which is in turn surrounded by a PU foot shell.

---

## 4 Maintenance

Maintenance must be carried out by competent personnel.

It is recommended that the following maintenance is carried out annually:

- Check for visual defects that may affect proper function
- Clean and relubricate bearing surface. Use Sapphire Endure (928015) or equivalent PTFE loaded grease (see Section 8 Maintenance Instructions).
- Remove the foot shell and glide sock, check for damage or wear and replace if necessary.
- Check all screws for tightness, clean and reassemble as necessary.
- Visually check the heel and toe springs for signs of delamination or wear and replace if necessary. Some surface damage may occur after a period of use, this does not affect the function or strength of the foot.

Any changes in performance of this device must be reported to the practitioner. Changes in performance may include:

- Increase in ankle stiffness
- Changes in axial or torsional stiffness
- Reduced ankle support (free movement)
- Softening of springs
- Excessive play in bearings
- Any unusual noise

The user should be advised to contact their practitioner if their condition changes.

The practitioner must also be informed of any changes in body weight and/or activity level.

The user should be advised that a regular visual check of the foot is recommended, signs of wear that may affect function should be reported to their service provider (e.g. significant wear or excessive discoloration from long term exposure to UV).

### **Cleaning:**

Use a damp cloth and mild soap to clean outside surfaces. Do not use aggressive cleansers.

---

## 5 Limitations on use

### **Intended life:**

A local risk assessment should be carried out based upon activity and usage.

### **Lifting loads:**

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

### **Environment:**

Avoid exposing EchelonVT to corrosive elements such as water, acids and other liquids. Also avoid abrasive environments such as those containing sand for example as these may promote premature wear.

Suitable for outdoor use


Exclusively for use between -15°C and 50°C (5°F to 122°F).

It is recommended that only Endolite products be used in conjunction with EchelonVT.


# 6 Bench alignment

## 6.1 Static Alignment

Align transfemoral devices according to fitting instructions supplied with the knee.  
Keep the build line between pivots as shown, using shift and/or tilt devices as necessary.


Cut shin to approximate length, allowing 10-15mm for final trimming. **Note:** trial length will typically include 10mm to allow for axial compression of heel, toe and spring (axial).

\* Allow for users own footwear

## 6.2 Biomimetic Alignment

The aim of alignment is to achieve a “balance point” while standing and set the hydraulically damped range of motion. The aim of damping adjustment is to fine tune the ankle-foot roll-over stiffness characteristics until a comfortable gait is achieved. Due to the increased range of motion provided by the ankle the user may experience the need for more voluntary control and initially find the ankle disconcerting during setup. This should quickly pass upon completion of satisfactory setup.


Ensure that the user is relaxed and not resting on the dorsi-flexion limit.

## 6.3 Biomimetic Adjustment

N.B. Carry out static alignment while ensuring the user has some means of support such as parallel bars. This is standing alignment only.


Use shift for static alignment and standing.

The device should encourage some degree of self adjustment to achieve a sense of balance for the user during standing.

## 6.4 Dynamic Adjustment

### Adjustment of the hydraulic valves

The user should experience the ankle moving with the body through the gait cycle. There should be no effort exerted by the user to overcome the hydraulic resistance of the ankle.

#### Guidance:

Following dynamic adjustment, trial the foot/ankle on ramps and stairs. Ensure the user is comfortable with the kind of terrain he/she may normally be expected to encounter. If the user reports any issues with comfort, usability or range of movement of the ankle, adjust accordingly.

## 7 Fitting Advice

The correct alignment (A-P position), range of motion (distribution of plantarflexion to dorsiflexion) and adjustment of the hydraulic settings are critical in achieving a smooth roll over and correct slope adaptation (see 6.3).

EchelonVT axial springs are designed to work together with the heel and toe springs to give smooth progression for most users. The axial spring controls both compressive and torsional resilience and typical axial movement should be between 3-6mm for normal walking.

EchelonVT is available with a choice of 4 axial spring rates (see section 1.0 for details).

Allow the user some time to become familiar with the initial set-up.

If after following the instructions below you still have problems with the function please contact the sales team in your area for advice.

Any of the following will have a negative effect on function and stability:

- Incorrect spring selection
- Incorrect A-P shift alignment
- Incorrect distribution of plantarflexion and dorsiflexion range

	<i>Symptoms</i>	<i>Remedy</i>
1.	<ul style="list-style-type: none"> <li>• Sinking at heel strike</li> <li>• Difficulty in achieving a smooth progression to mid stance</li> <li>• User feels they are walking up hill or forefoot feels excessively long</li> </ul>	<ul style="list-style-type: none"> <li>• Increase plantarflexion resistance</li> <li>• Check A-P shift alignment; ensure foot is not too anteriorly positioned</li> <li>• Check distribution of plantarflexion and dorsiflexion movement; ensure that the plantarflexion range is not excessive</li> <li>• Check spring category is not too soft, if so fit higher rate spring</li> </ul>
2.	<ul style="list-style-type: none"> <li>• Progression from heel strike to mid stance is too rapid</li> <li>• Difficulty in controlling the energy return from the foot at the heel strike (reduced knee stability)</li> <li>• User feels heel is too hard, fore foot is too short</li> </ul>	<ul style="list-style-type: none"> <li>• Reduce plantarflexion resistance</li> <li>• Check A-P shift alignment; ensure foot is not too posteriorly positioned</li> <li>• Check distribution of plantarflexion and dorsiflexion movement; ensure that there is adequate plantarflexion range</li> <li>• Check the spring category is not too high for the weight and activity of the patient, if so fit lower rate spring</li> </ul>

	<i>Symptoms</i>	<i>Remedy</i>
3.	<p>Heel contact and progression feel OK but:</p> <ul style="list-style-type: none"> <li>• Forefoot feels too soft</li> <li>• Forefoot feels too short</li> <li>• User feels they are walking <i>down</i> hill, possibly with reduced knee stability</li> <li>• Lack of energy return</li> </ul>	<ul style="list-style-type: none"> <li>• Increase dorsiflexion resistance</li> <li>• Check A-P shift alignment; ensure foot is not too posteriorly positioned</li> <li>• Check distribution of plantarflexion and dorsiflexion movement; ensure that there is not excessive dorsiflexion range</li> <li>• Check the spring category is not too soft for the weight and activity of the patient, if so fit higher rate spring</li> </ul>
4.	<p>Heel contact and progression feel OK but:</p> <ul style="list-style-type: none"> <li>• Forefoot feels too rigid</li> <li>• Forefoot feels too long</li> <li>• Feels like walking up hill</li> </ul>	<ul style="list-style-type: none"> <li>• Reduce dorsiflexion resistance</li> <li>• Check A-P shift alignment; ensure foot is not too anteriorly positioned</li> <li>• Check distribution of plantarflexion and dorsiflexion movement; ensure that there is sufficient dorsiflexion range</li> <li>• Check the spring category is not too rigid for the weight and activity of the patient, if so fit lower rate spring</li> </ul>

# 8 Maintenance Instructions

**Be aware of finger trap hazard at all times.**

## Foot Shell Removal

**1**

Insert shoe horn behind heel spring.

**2**

Rotate shoe horn as shown to remove shell.

## Heel & Toe Spring Replacement


**3**

Remove heel spring and screws.

**4**

Remove toe spring screw, replace toe. Upon reassembly, use Loctite 243 (926012) and torque to 35Nm. Ensure toe spring is central to the carrier.

**5**

Re-assemble with replacement heel spring. Use Loctite 243 (926012) and torque to 15Nm.

**6**

Cover appropriate lines on carrier with permanent black marker to leave springset number showing.

## 8 Maintenance Instructions *(cont.)*

7

If a foam cosmesis is to be fitted, roughen top surface of foot shell to provide ideal bonding surface.

8

Lubricate toe and heel if required.  
(Foot shell is pre-lubricated).


9

toe spring location in foot shell.


10

Use a suitable lever to encourage the heel spring into location in the foot shell.

11



Ensure heel spring is engaged into slot.

If a cosmetic finish is required please contact a member of the Endolite Sales Team.

## 8 Maintenance Instructions *(cont.)*

### Pylon Sleeve

Clean and relubricate bearing surface.

Use Sapphire Endure (928015) or equivalent PTFE loaded grease.

**Re-assembly:**

Ensure the rotation stop screw aligns with the slot in the pylon apply Loctite 243 and tighten.


Replacement Gaiter available Part No. 534003

# 9 Technical Data

Operating and Storage Temperature Range:	-15°C to 50°C (5°F to 122°F)
Component Weight (Size 26N):	1.05kg ( 2lb 3ozs)
Recommended Activity Level:	3
Maximum User Weight:	125kg (275lb)
Sizes:	22 to 30cm
Proximal Alignment attachment:	Male Pyramid (Endolite) +/- 7°
Range of Ankle Movement:	6° plantar to 3° dorsi flexion
Maximum Rotation:	+/- 15°
Typical Axial Movement:	3- 6mm [maximum 10mm]
Build Height:*	22cm - 24cm 168mm 25cm - 26cm 173mm 27cm - 30cm 178mm
Heel Height:	10mm

## Fitting Length

**\* Note:**  
Refer to  
Section 6

Size	A
22-24	168 mm
25-26	173 mm
27-30	178 mm

Size	B
22-26	65 mm
27-28	70 mm
29-30	75 mm

# 10 Ordering Information

Spring Kits				
Rate	Foot Sizes			
	Small (S)	Medium (M)	Large (L)	Extra Large (XL)
	22-24	25-26	27-28	29-30
Set 1	539801S	539810S	539819S	539828S
Set 2	539802S	539811S	539820S	539829S
Set 3	539803S	539812S	539821S	539830S
Set 4	539804S	539813S	539822S	539831S
Set 5	539805S	539814S	539823S	539832S
Set 6	539806S	539815S	539824S	539833S
Set 7	539807S	539816S	539825S	539834S
Set 8	539808S	539817S	539826S	539835S

Foot Shell (for dark add 'D')		
Size/Side	Narrow	Wide
22L	539038S	-
22R	539039S	-
23L	539040S	-
23R	539041S	-
24L	539042S	-
24R	539043S	-
25L	539044SN	539044SW
25R	539045SN	539045SW
26L	539046SN	539046SW
26R	539047SN	539047SW
27L	539048SN	539048SW
27R	539049SN	539049SW
28L	-	539050S
28R	-	539051S
29L	-	539052S
29R	-	539053S
30L	-	539054S
30R	-	539055S

Item	Part. No.
Glide Sock (Sizes 22-26)	531011
Glide Sock (Sizes 27-30)	532811
DF/PF Adjuster Key, 4 A/F Allen	940236
Gaiter	534003
Rotation stop screw	534005

**Liability**

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for damage caused by the component combinations that were not authorized by the manufacturer.

**CE Conformity**

This product meets the requirements of 93/42/EEC guidelines for medical products. This product has been classified as a Class 1 Product according to the classification criteria outlined in Appendix IX of the guidelines. The Declaration of Conformity was therefore created by Blatchford Products Limited with sole responsibility according to Appendix VII of the guidelines.

**Warranty**

This device is warranted for 36 months - foot shell 12 months - glide sock 3 months. The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions. See Endolite website for the current full warranty statement.

**Environmental Aspects**

Where possible the components should be recycled in accordance with local waste handling regulations.

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Endolite India Ltd
Endolite Deutschland GmbH,